ECOMED® ANTIMICROBIAL & ANTIVIRAL DISPOSABLE HOSPITAL CURTAINS

SURROUND YOU AND PROTECT YOU.

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TEO HONG SILOM CO., LTD.

Tel : 02 312 0045-69 # 1709, 1720, 1743

Email : med_mkt@teohong.com

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*Testing reports are available upon request.

Bacteria	International Standards	Log Reduction
E. coli	ISO 20743	4.40
S. aureus	ISO 20743	3.60
MRSA	ISO 20743	4.40
VRE	ISO 20743	3.30
Ps. aeruginosa	ISO 20743	4.30
E. faecalis	ISO 20743	4.40
A. baumannii	ISO 20743	4.50
K. pneumoniae	ISO 20743	5.90
S. enteritidis	ISO 20743	5.60

Viruses	Antiviral Active Value
Human Coronavirus (HCoV-229E)	99.99%
H1N1 Influenza A virus (Swine Flu)	99.99%
Human enterovirus 71 (EV-A71)	99.99%
Coxsackievirus A21	99.63%
Adenovirus type 3	99.56%

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Country	International Standards
United Kingdom	BS 5867 Part 2 Type C: 2008
United States	NFPA 701: 2023
Europe	EN13772, EN13773 Class 1
Australia	AS 2755.2-1985
New Zealand	AS/NZS 1530.3: 1999

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QUALITY ASSURANCE PROGRAM

QA TEST CERTIFICATION

This certifies that the products listed below have been tested and testing has been deemed successful per the requirements of the test method.

Changzhou EcoMed Technologies Co., Ltd.

IMSL Certificate #1046706.1A, 1B, 1C

TEST METHOD	SAMPLE DESCRIPTION	TEST DATE	EFFICACY ESTABLISHED
ISO 20743	PP Non-woven medical curtain with 0.5% CA-16	29 NOVEMBER 2023	Antibacterial: E. coli S. aureus Ps. aeruginosa MRSA VRE

Veronica Steinmetz Digitally signed by

Steinmetz Date: 2023.12.19

SIGNED: DATE: DECEMBER 19, 2023

Veronica Steinmetz

Quality Assurance Specialist

MICROBAN PRODUCTS COMPANY

11400 Vanstory Drive | Huntersville | NC 28078 | United States | T: +1 (704) 875-0806

MICROBAN PRODUCTS COMPANY

QUALITY ASSURANCE PROGRAM

QA TEST CERTIFICATION

This certifies that the products listed below have been tested and testing has been deemed successful per the requirements of the test method.

Changzhou EcoMed Technologies Co., Ltd.

IMSL Certificate #1047460.1A-1B-1C-1D

TEST METHOD	SAMPLE DESCRIPTION	TEST DATE	EFFICACY ESTABLISHED
MOD ISO 20743	PP Non-woven Hospital Curtain (0.5% CA-16): 3 year used curtain - Batch E01-Apr2021	28 June 2024	Antibacterial: E.coli S.aureus Ps. aeruginosa MRSA VRE A.baumannii

Veronica Steinmetz Digitally signed by
Veronica Steinmetz
Date: 2024.07.17
10:20:30 -04'00'
SIGNED: _____ DATE: JULY 17, 2024

Veronica Steinmetz
Quality Assurance Specialist

MICROBAN PRODUCTS COMPANY
11400 Vanstory Drive | Huntersville | NC 28078 | United States | T: +1 (704) 875-0806

CERTIFICATE OF ANALYSIS**METHOD: (Determination of Antibacterial Activity using Modified ISO 20743:2021)****CUSTOMER:**

Microban (Europe) Ltd
Pera Business Park
Nottingham Road
Melton Mowbray
Leicestershire
LE13 1HA

Certificate No: 1046706.1A**Date Received:** 21-Nov-23**Date Analysed:** 29-Nov-23**Date Reported:** 01-Dec-23**MBE Ref:** QA-2023**SAMPLE DETAILS****Customer Ref:** Changzhou Ecomed Technologies Co Ltd**DEVIATIONS FROM STANDARD**

The inoculum was prepared from plate culture in 1: 500 Nutrient Broth as described in ISO 22196 : 2011

Escherichia coli was used in place of *Klebsiella pneumoniae* as the gram negative species

Each replicate was prepared as 0.4 g test pieces and inoculated with 0.2 ml of the relevant inoculum, with the addition of a 0.05% wetting agent. Efficacy is calculated as CFU cm⁻²

TEST RESULTS (AS CFU CM⁻²)

SAMPLE	SPECIES	CONTACT TIME (HRS)		LOG REDUCTION	% DIFFERENCE
		0	24		
IMSL Control Fabric	<i>E. coli</i>	2.2E+05	3.6E+07	-	-
PP NW Hospital Curtain (0.5% CA-16)	<i>E. coli</i>	2.2E+05	5.4E+01	3.6	99.98%
IMSL Control Fabric	<i>S. aureus</i>	2.7E+05	2.2E+05	-	-
PP NW Hospital Curtain (0.5% CA-16)	<i>S. aureus</i>	2.7E+05	8.2E+02	2.5	99.70%

The above data shows the difference in the population following contact with the surface of the samples listed for 24 hours at 35°C under a Relative Humidity (RH) of > 95% relative to the initial population.

INDUSTRIAL MICROBIOLOGICAL SERVICE LTD
PALE LANE
HARTLEY WINTNEY
HANTS
RG27 8DH
UK

RICHARD WEBB
SCIENTIST

Industrial Microbiological Services Ltd Registered in England No 3284623 Registered Office The Oldfellow Hall Oxford Road Reading Berkshire RG1 2NG

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CERTIFICATE OF ANALYSIS**METHOD:** (Determination of Antibacterial Activity using Modified ISO 20743:2021)**CUSTOMER:**

Microban (Europe) Ltd
Pera Business Park
Nottingham Road
Melton Mowbray
Leicestershire
LE13 1HA

Certificate No: 1046706.1B**Date Received:** 21-Nov-23**Date Analysed:** 29-Nov-23**Date Reported:** 01-Dec-23**MBE Ref:** QA-2023**SAMPLE DETAILS****Customer Ref:** Changzhou Ecomed Technologies Co Ltd**DEVIATIONS FROM STANDARD**

The inoculum was prepared from plate culture in 1: 500 Nutrient Broth as described in ISO 22196 : 2011
Pseudomonas aeruginosa was used in place of *Klebsiella pneumoniae* as the gram negative species and
Methicillin Resistant *Staphylococcus aureus* in place of *Staphylococcus aureus* as the gram positive species
Each replicate was prepared as 0.4 g test pieces and inoculated with 0.2 ml of the relevant inoculum, with the
addition of a 0.05% wetting agent. Efficacy is calculated as CFU cm⁻²


TEST RESULTS (AS CFU CM⁻²)

SAMPLE	SPECIES	CONTACT TIME (HRS)		LOG REDUCTION	% DIFFERENCE
		0	24		
IMSL Control Fabric	<i>Ps aeruginosa</i>	2.4E+05	6.5E+07	-	-
PP NW Hospital Curtain (0.5% CA-18)	<i>Ps aeruginosa</i>	2.4E+05	≤ 12.5	≥ 4.3	≥ 99.99%
IMSL Control Fabric	<i>MRSA</i>	2.1E+05	6.6E+04	-	-
PP NW Hospital Curtain (0.5% CA-18)	<i>MRSA</i>	2.1E+05	4.9E+01	3.6	99.98%

The above data shows the difference in the population following contact with the surface of the samples listed for 24 hours at 35°C under a Relative Humidity (RH) of > 95% relative to the initial population.

INDUSTRIAL MICROBIOLOGICAL SERVICE LTD
PALE LANE
HARTLEY WINTNEY
HANTS
RG27 8DH
UK

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CERTIFICATE OF ANALYSIS**METHOD: (Determination of Antibacterial Activity using Modified ISO 20743:2021)****CUSTOMER:**

Microban (Europe) Ltd
Pera Business Park
Nottingham Road
Melton Mowbray
Leicestershire
LE13 1HA

Certificate No:

1046706.1C

Date Received:

21-Nov-23

Date Analysed:

29-Nov-23

Date Reported:

01-Dec-23

MBE Ref:

QA-2023

SAMPLE DETAILS**Customer Ref:**

Changzhou Ecomed Technologies Co Ltd

DEVIATIONS FROM STANDARD

The inoculum was prepared from plate culture in 1: 500 Nutrient Broth as described in ISO 22196 : 2011

Vancomycin-resistant Enterococcus was used in place of *Staphylococcus aureus* as the gram negative species

Each replicate was prepared as 0.4 g test pieces and inoculated with 0.2 ml of the relevant inoculum, with the addition of a 0.05% wetting agent. Efficacy is calculated as CFU cm⁻²

TEST RESULTS (AS CFU CM⁻²)

SAMPLE	SPECIES	CONTACT TIME (HRS)		LOG REDUCTION	% DIFFERENCE
		0	24		
IMSL Control Fabric	VRE	2.6E+05	1.9E+05	-	-
PP NW Hospital Curtain (0.5% CA-16)	VRE	2.6E+05	≤ 12.5	≥ 4.3	≥ 99.99%

The above data shows the difference in the population following contact with the surface of the samples listed for 24 hours at 35°C under a Relative Humidity (RH) of > 95% relative to the initial population.

INDUSTRIAL MICROBIOLOGICAL SERVICE LTD
PALE LANE
HARTLEY WINTNEY
HANTS
RG27 8DH
UK

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Date: 2024-08-26

Sample Description:

Sample No.	SGS Sample ID	Description
1	GZF24-0034309-0001	sample in bag

Test Result(s):

Test Requested: Selected test(s) as requested by applicant: Antimicrobial activity test

Test Method: Refer to ISO 20743:2021(E) Textiles-Determination of antibacterial activity of textile products: 8.1 Absorption method

GZF24-0034309-0001

Test organism(s)	Klebsiella pneumoniae ATCC 4352
Test inoculum (CFU/mL)	1.8×10^5
F	2.6
lgC _i	7.2
lgC ₀	4.6
G	-3.2
lgT _i	1.3
lgT ₀	4.5
A*	5.9
Antibacterial activity(%)***	>99

Notes:

The control sample was 100% cotton fabric, provided by SGS lab.

F: Growth value of F on the control specimens ($F = \lg C_i - \lg C_0$);

G: Growth value of G on the antibacterial testing specimens ($G = \lg T_i - \lg T_0$);

A*: antibacterial activity value ($A = F - G$);

lgC_i: common logarithm of arithmetic average of the numbers of bacteria, obtained from three control specimens after 24 h incubation;

lgC₀: common logarithm of arithmetic average of the numbers of bacteria, obtained from three control specimens immediately after incubation;

lgT_i: common logarithm of arithmetic average of the numbers of bacteria, obtained from three test specimens after 24 h incubation;

lgT₀: common logarithm of arithmetic average of the numbers of bacteria, obtained from three test specimens immediately after incubation.

The testing specimens sterilization method was UV sterilization of both sides for 30 min.

The testing specimens incubation time was 24 h.

The testing measuring method was plate count method.

* In case of $C_0 > T_0$, substitute C_0 for T_0 .

*** The percentage of Antibacterial activity value(%) is added by client requirement, the calculate formulation is $[(C_i - T_i) / C_i] \times 100$.

According to client's requirement, the inoculum was prepared from plate culture in 1: 500 Nutrient Broth.

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Date: 2024-05-28

Sample Description:

Sample No.	SGS Sample ID	Description
1	ASH24-0027972-0001	sample in bag

Test Result(s):

Test Requested: Selected test(s) as requested by applicant

Test Method: ISO 18184:2019 Textiles-Determination of antiviral activity of textile products

ASH24-0027972-0001

Test virus strain Virus strain and host cell Influenza A virus H1N1 ATCC-VR-1469, MDCK cell ATCC-CCL-34

	Logarithm of infectivity titre value immediate after inoculation of the control specimen(IgTCID ₅₀ / vial)	Logarithm of infectivity titre value after 2h contacting with the control specimen(IgTCID ₅₀ / vial)	Logarithm of infectivity titre value after 2h contacting with the antiviral fabric specimen(IgTCID ₅₀ / vial)
sample 1	6.68	6.30	2.80
sample 2	6.55	6.43	2.80
sample 3	6.68	6.55	2.80
Average logarithm of infectivity titre value (IgTCID ₅₀ /vial)	6.64	6.43	2.80
Average infectivity titre value (TCID ₅₀ /vial)	4.37×10^6	2.74×10^6	6.31×10^2
Mv:Antiviral activity value		3.84	
Antiviral activity rate(%) [*]		99.99	

Notes:

- 1.The method to determine the infectious titre: TCID₅₀ method.
 - 2.Control fabric: cotton 100% woven fabric provided by lab.
 - 3.Sterilization of antiviral fabric specimen: UV 30min
 - 4.Contact time of the inoculated specimens:2h
 - 5.The concentration of the virus inoculum suspension: 3.16×10^7 TCID₅₀/mL
- ^{*} Antiviral activity rate(%) calculate formulation : $[(A-C) / A] * 100\%$
A: Average infectivity titre value immediate after inoculation of the control specimen(TCID₅₀/ vial)
C: Average infectivity titre value after 2h contacting with the antiviral fabric specimen(TCID₅₀/ vial)

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MBE Ref: QA 2024

SAMPLE DETAILS

Customer Ref: Changzhou EcoMed Technologies Co., Ltd.

DEVIATIONS FROM STANDARD

The inoculum was prepared from plate culture in 1: 500 Nutrient Broth as described in ISO 22196 : 2011
Escherichia coli, *Pseudomonas aeruginosa* and *Acinetobacter baumannii* was used in place of *Klebsiella pneumoniae* as the gram negative species. Each replicate was prepared as 0.4 g test pieces and inoculated with 0.2 ml of the relevant inoculum, with the addition of a 0.05% wetting agent. Efficacy is calculated as CFU g⁻¹

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INDUSTRIAL MICROBIOLOGICAL SERVICE LTD
PALE LANE
HARTLEY WINTNEY

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MBE Ref: QA 2024

SAMPLE DETAILS

Customer Ref: Changzhou EcoMed Technologies Co., Ltd.

DEVIATIONS FROM STANDARD

The inoculum was prepared from plate culture in 1: 500 Nutrient Broth as described in ISO 22196 : 2011
Escherichia coli, *Pseudomonas aeruginosa* and *Acinetobacter baumannii* was used in place of *Klebsiella pneumoniae*
as the gram negative species. Each replicate was prepared as 0.4 g test pieces and inoculated with 0.2 ml of the
relevant inoculum, with the addition of a 0.05% wetting agent. Efficacy is calculated as CFU g⁻¹

TEST RESULTS (AS CFU G⁻¹)

SAMPLE	SPECIES	CONTACT TIME (HRS)		LOG REDUCTION	% DIFFERENCE
		0	24		
IMSL Control Fabric	<i>E. coli</i>	3.1E+05	1.7E+07	-	-
Non Woven Curtain - Control	<i>E. coli</i>	3.1E+05	3.5E+06	-	-
Non-woven Curtain : 0.5% CA-16 : Batch E01-Apr2021 : 3 Years	<i>E. coli</i>	3.1E+05	≤ 12.5	≥ 4.4	≥ 99.99%
IMSL Control Fabric	<i>S. aureus</i>	3.2E+05	5.4E+05	-	-
Non Woven Curtain - Control	<i>S. aureus</i>	3.2E+05	6.0E+03	-	-
Non-woven Curtain : 0.5% CA-16 : Batch E01-Apr2021 : 3 Years	<i>S. aureus</i>	3.2E+05	9.0E+01	3.6	99.97%
IMSL Control Fabric	VRE	2.0E+05	1.8E+05	-	-
Non Woven Curtain - Control	VRE	2.0E+05	2.9E+04	-	-
Non-woven Curtain : 0.5% CA-16 : Batch E01 Apr2021 : 3 Years	VRE	2.0E+05	1.1E+02	3.3	99.95%

The above data shows the difference in the population following contact with the surface of the samples listed
for 24 hours at 35°C under a Relative Humidity (RH) of > 95% relative to the initial population.

INDUSTRIAL MICROBIOLOGICAL SERVICE LTD
PALE LANE
HARTLEY WINTNEY

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EU Declaration of Conformity NON-STERILIZED DISPOSABLE MEDICAL CURTAINS

Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 concerning Medical Devices. The undersigned declares, under their sole responsibility, that the products described in this document meet the Council provisions that apply to them and the CE Mark may be affixed.

General Product Name	NON-STERILIZED DISPOSABLE MEDICAL CURTAINS
Legal Manufacturer:	Changzhou EcoMed Technologies Co., Ltd. Business Operation Address: No. 2 Jianye Road, New District, Changzhou City, Jiangsu Province, P.R. China
Description:	Single use medical curtains packed as bulk non-sterile units.
Intended Purpose:	Disposable medical devices for infection control and prevention.
MDR Classification:	Class 1 as per Rule 1 of Annex VIII
EU Authorised Representative:	Carotechion GmbH Niederrheinstraße. 71, 40474 Düsseldorf, Germany
EU Authorised Representative SRN:	DE-AR-000005946
Medical Device Regulation Assessment Route:	Self-certification by Medical Device Regulation Article 52: Conformity assessment based on technical documentation (Annex II & Annex III), and EU Declaration of Conformity (Annex IV).

Applicable Standards

This present declaration is also in conformity with the following European standards and Common Specifications (CS): 2017/745, ISO 13485:2016+A11:2021, ISO 14971:2019+A11:2021, ISO 15223-1:2021, ISO 20417:2021, BS EN ISO 10993 -1:2020.

Declaration of Conformity

Changzhou EcoMed Technologies Co., Ltd., with sole responsibility for the manufacture, packaging and labelling of the NON-STERILIZED DISPOSABLE MEDICAL CURTAINS listed on this certificate, declares that they conform to the relevant provisions of Regulation (EU) 2017/745 and ISO 13485 registered Quality Management System, regardless of whether these operations are carried out by the manufacturer or on their behalf by a third party.

Date: 18 JANUARY, 2024 Place: CHANGZHOU, CHINA Position: CEO Signature:

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